

X-RAY PROTECTIVE CLOTHING

LytaType jacket (vest) model 915

Velcro fastened overlapping front provides all round protection.
Front doubles the stated protection at an overlap of approximately 15 cm.

| Lead equivalent Front and back | Length | Style | Product code |
|---|--------|------------|--------------|
| Small for chest sizes 80 – 90 cm | | | |
| 0.25 / 0.25 mm | 52 cm | Slim fit | 915/22/S |
| 0.25 / 0.35 mm | 52 cm | " | 915/23/S |
| 0.35 / 0.35 mm | 52 cm | " | 915/33/S |
| Medium for chest sizes 90 – 104 cm | | | |
| 0.25 / 0.25 mm | 58 cm | Medium fit | 915/22/M |
| 0.25 / 0.35 mm | 58 cm | " | 915/23/M |
| 0.35 / 0.35 mm | 58 cm | " | 915/33/M |
| Large for chest sizes 104 – 120 cm | | | |
| 0.25 / 0.25 mm | 65 cm | Large fit | 915/22/L |
| 0.25 / 0.35 mm | 65 cm | " | 915/23/L |
| 0.35 / 0.35 mm | 65 cm | " | 915/33/L |

Optional extras

Pocket without embroidery

Embroidery in a position to be specified

900/001

9/B20



Models 915 & 919
Pockets 900/001 shown
are **optional extras**

LytaType colours and codes



Purple
01



Royal
02



Navy
04



Poppy
11



Burgundy
12



Silver grey
40



Bottle green
41



915

Please note the colours shown are approximate. A colour swatch is available on request.

Tailor made *LytaType* personal x-ray protective clothing

Please contact us for required measurements

ORDER GUIDE. Please specify: -

- | | | | |
|---|-----------------|-----------------------|---------|
| 1. Model description | 2. Model number | 3. Lead equivalent | 4. Size |
| 5. Length, taken from the middle shoulder | 6. Colour | 7. Any optional extra | |