

X-RAY PROTECTIVE CLOTHING

LytaType wrap-around coat model 905

Front opening all-round coat with front Velcro fastening, fitted with padded shoulders for extra comfort.
Wide front doubles the stated protection at the overlap.

Lead equivalent Front and back

Length

Product code

Small for chest sizes 80 – 90 cm

Protection is double where the front of the apron overlaps.

Lead equivalent	Length	Fit	Product code
0.25 / 0.25 mm	85 cm	Slim fit	905/22/085/S
"	90 cm	"	905/22/090/S
"	100 cm	"	905/22/100/S
"	110 cm	"	905/22/110/S
0.25 / 0.35 mm	85 cm	"	905/23/085/S
"	90 cm	"	905/23/090/S
"	100 cm	"	905/23/100/S
"	110 cm	"	905/23/110/S
0.35 / 0.35 mm	85 cm	"	905/33/085/S
"	90 cm	"	905/33/090/S
"	100 cm	"	905/33/100/S
"	110 cm	"	905/33/110/S

Medium for chest sizes 90 – 104 cm

Protection is double where the front of the apron overlaps.

Lead equivalent	Length	Fit	Product code
0.25 / 0.25 mm	90 cm	Medium fit	905/22/090/M
"	100 cm	"	905/22/100/M
"	110 cm	"	905/22/110/M
0.25 / 0.35 mm	90 cm	"	905/23/090/M
"	100 cm	"	905/23/100/M
"	110 cm	"	905/23/110/M
0.35 / 0.35 mm	90 cm	"	905/33/090/M
"	100 cm	"	905/33/100/M
"	110 cm	"	905/33/110/M

Large for chest sizes 104 – 120 cm

Protection is double where the front of the apron overlaps.

Lead equivalent	Length	Fit	Product code
0.25 / 0.25 mm	90 cm	Large fit	905/22/090/L
"	100 cm	"	905/22/100/L
"	110cm	"	905/22/110/L
0.25 / 0.35 mm	90 cm	"	905/23/090/L
"	100 cm	"	905/23/100/L
"	110 cm	"	905/23/110/L
0.35 / 0.35 mm	90 cm	"	905/33/090/L
"	100 cm	"	905/33/100/L
"	110 cm	"	905/33/110/L

Also available: -

Front 0.25 mm, back 0.50 mm @ 0.25 / 0.25 mm prices <u>plus</u> 20%	905/25/.../.
Front 0.35 mm, back 0.50 mm @ 0.25 / 0.25 mm prices <u>plus</u> 25%	905/35/.../.
Front 0.50 mm, back 0.50 mm @ 0.25 / 0.25 mm prices <u>plus</u> 30%	905/55/.../.
Please add length and S, M, or L (small, medium or large) to the applicable product code.	
Example: 0.25 front (2), 0.50 back (5), 90 cm long and medium size:	905/25/090/M

LytaType wrap-around coat model 903

Front opening all-round coat with front Velcro fastening, fitted with padded shoulders for extra comfort.

Narrow front doubles the stated protection at an overlap of approximately 15 cm.

Sizing is the same as for model 905

Optional extras

Pocket	900/001
Embroidery in a position to be specified	9/B20
Lumbar support belt elasticated type transfers weight to the hips	9/B15

Please contact us for required measurements.



Model 905 with a wide front overlap

Optional extras shown: -

- Pocket 900/001
- 9/B15 support belt,
- 951/05 thyroid collar

ORDER GUIDE Please specify:

- Model description
- Model number
- Lead equivalent
- Size
- Length taken from the middle shoulder
- Colour (see page 13)
- Any optional item

LytaType colours and codes



Purple
01



Royal
02



Navy
04



Poppy
11



Burgundy
12



Silver grey
40



Bottle green
41

Please note the colours shown are approximate. A colour swatch is available on request.

Tailor made *LytaType* personal x-ray protective clothing

Please contact us for required measurements.