

Atomscope HF400A

The high frequency, high power 6kW Atomscope HF400A moves portable x-ray units to a powerful new level

- > KV range from 40 to 120 in 2kV steps
- > Wipe clean touch panel displays
- > Digital display of kV, mAs and time
- > LED display "X-RAY", "READY", "ERROR"
- > Hand exposure switch for collimator light
- > Acoustic and optical timer signal
- > Tape measure for SID
- > Thickness gauge for measuring object thickness
- > Switchable audible exposure warning
- > 5 preset memory functions
- > Powered by 13 amp domestic supply
- > Aluminum casing
- > 5 year extended warranty



Euromed Medical Systems

Ph 09 411 7508

Email: thomas.heus@euromed.co.nz